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i3 FOUR - inline sterile filter for secure germ retention

Intended use


The i3 FOUR inline filter is used for the sterile filtration of tap water to prevent contamination due to waterborne bacteria such as *Legionella pneumophila*, *Pseudomonas aeruginosa* and mycobacteria. It is installed in the inflow pipes of devices and instruments. The areas of use of the i3 FOUR are varied and include many different situations. For example, it is used in medical settings for the reprocessing of devices (e.g. endoscopes) and in dental and ENT units. It provides reliable protection against waterborne pathogens in the provision of drinking water (drinking water fountains), air conditioning, laboratory equipment and other situations.

Application

The i3 FOUR inline sterile filter is mounted in the pipe system by means of quick-release connectors. After a 50-day maximum period of use the disposable filter must be dismantled and disposed of. Please read the operating instructions thoroughly before use.

Safety instructions

- Installation of the filter and adapters may only be carried out by qualified technical personnel
- The water supply must be closed when the filter is not in use
- There must not be any tension (leverage of the connection hoses) on the installed filter
- Do not install the filter/accessories in tankless (low pressure) systems
- Do not use the filter outside of its specification
- The filter can be used for a maximum of 50 days
- The duration of use may be shorter if water quality is poor (rust, sand, etc.)
- It is not intended to be used for injection, infusion, irrigation or hemodialysis
- The filter does not change the chemical properties nor the taste or smell of water
- Do not remove the filter while it is in use – air in the piping system can block the filter
- If there are any signs of contamination, then the water filter should be replaced

 Article number	 Latex free	 Do not use if the filter or the packaging is damaged
 Read the instructions thoroughly before use	 Serial number	
 Do not use after the expiry date	 LOT number	 Manufacturer

Preparation

When installing the filter for the first time, a suitable quick release coupling with seal must be fitted (see i3 accessories). Do not use the filter if the packaging is damaged. Do not use the filter if the use-by date has elapsed (see filter packaging).

Installation

1. Disinfect your hands
2. Perform a surface disinfection of the quick release coupling
3. Open the packaging of the filter. If possible, do not touch the filter inlet and outlet side
4. If possible, close the water supply and insert the outlet side of the filter into the corresponding quick release coupling (see arrow on product).
5. Install the inlet side of the filter on the second quick release coupling
6. Check whether the filter is firmly installed in the quick release couplings.
Open the water supply and check the installation for leaks

Deinstallation

7. If possible, close the water supply and disconnect the filter on the inlet side
8. Then disconnect the filter on the outlet side
9. Dispose the filter with household trash or according to the local guidelines

Technical data

Bacterial retention:	> 7 log (<i>Brev. dim.</i>) per cm ² filter area
Flow rate:	6 l/min at 5 bar
Duration of use:	50 days (7 weeks)
Maximum operating pressure:	5 bar
Maximum operating temperature:	60 °C (max. 70 °C for 30 mins.)
Storage life:	5 years (at max. 40° C storage temperature)
Chlorine resistance:	1 hour ≤ 200 ppm active chlorine; 50 days ≤ 10 ppm active chlorine

Disinfection:	Can be disinfected by wiping with common surface disinfectants 41 x 82 mm (diameter x length) Quick release coupling
Measurements:	
Adaption:	Quick release coupling

Accessories

G3/8" AG Inline-Quick release coupling (Art. No. 2060005)
(additional accessories on request)

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G3/8" AG Inline-Snelkoppeling (Art. Nr. 2060005)
i3 Membrane GmbH